

E.Z.N.A.[®] FFPE & Blood

	<u>Blood RNA Kit</u> <u>R6814</u>	<u>PX Blood RNA Kit</u> <u>R1057</u>	<u>Blood RNA Midi Kit</u> <u>R6615</u>	<u>Blood RNA Maxi Kit</u> <u>R6616</u>	<u>Mag Bind Blood RNA Kit</u> <u>M2837</u>	<u>FFPE RNA Kit</u> <u>R6954</u>	<u>Mag Bind FFPE RNA Kit for KF</u> <u>M6953</u>	<u>E-Z 96 Mag Bind FFPE RNA Kit</u> <u>M2551</u>
Starting Sample Volume	10-1000 µl	1 Paxgene™ Tube	2 mL-10 mL	10 mL-50 mL	50 µl	3-8 Sections	3-8 Sections	3-8 Sections
Binding Format	Mini Spin Columns	Mini Spin Columns	Midi Spin Columns	Maxi Spin Columns	Magnetic Beads	Mini Spin Columns	Magnetic Beads	Magnetic Beads
Homogenizer Columns Included	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>				
Red Blood Cell Lysis Step required	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>				
Includes DNase Treatment		<input checked="" type="checkbox"/>			<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		
Organic Extractions Required						<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Automatable					<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>