

E.Z.N.A. Plasmid Mega and Giga

	<u>Fastfilter</u> <u>Plasmid Mega</u> D6929	<u>Endo-Free</u> <u>Plasmid Mega</u> D6228	<u>Plasmid</u> <u>Giga</u> D6920	<u>Endo-Free</u> <u>Plasmid Giga</u> D6234
Processing Time (Minutes)	40 min	60 min	60 min	80 min
Plasmid Type	High & Low Copy Plasmids	High & Low Copy Plasmids	High & Low Copy Plasmids	High & Low Copy Plasmids
Starting Amounts	200-500 mL	1000-2500 mL	200-500 mL	1000-2500 mL
Lysate Clearance Syringe	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		
Vacuum Manifold Required	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Restriction Enzyme, Sequencing	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Transfection and PCR	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Sensitive Cell Line & Primary Cell Transfection		<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>